



SHIJIAZHUANG HONGRAY GROUP

South Tongda Rd., East Dist. Jinzhou City Hebei, 052260, China

PRODUCT SPECIFICATION SHEET

GLOVE SPECIFICATION

Product Specification Group Code: NPF2001-2005, XS-XL

Product Name: Nitrile Examination Glove, Powder Free, Chlorinated, Non-sterile, Medical Grade/PPE Grade, M/4.5g

GENERAL INFORMATION

Type: Single-use disposable examination or protective glove

Shape: Ambidextrous

Material: Synthetic Nitrile Latex

Cuff: Rolled cuff

Surface: Full Textured

Thickness:

Single wall according to ASTM D6319, ASTM D3767	Palm: Min. 0.08mm	Finer: Min. 0.10mm
---	-------------------	--------------------

Residual Powder: (according to EN 455-3)	Max. 2.0mg/glove
---	------------------

Force at break: (according to EN 455-2)	Before Aging: Med. $\geq 6.0N$ After Aging: Med. $\geq 6.0N$
--	---

Length: (according to EN 455-2)	XS-XL: Med. $\geq 240mm$
------------------------------------	--------------------------

Tensile Strength: (according to ASTM D 6319, ASTM D412, ASTM D573)	Before Aging: Min. 14Mpa After Aging: Min. 14Mpa
--	---

Elongation %: (according to ASTM D 6319, ASTM D412, ASTM D573)	Before Aging: Min. 500% After Aging: Min. 400%
--	---

Available Sizes & Palm width (mm): (according to EN 455-2)

XS: 70 \pm 10 S: 80 \pm 10 M: 95 \pm 10 L: 110 \pm 10 XL: 120 \pm 10
--

Storage: 3 years storage life (date of expiration see on the box). Store in a cool, dry place. Shield open box from direct sunlight, fluorescent lighting or x-rays.

DOCUMENT NO.: SETINO-NG-002

INITIALED DATE: 2021.10.25

REVISION DATE: N/A

REVISION TIMES: 0



SHIJIAZHUANG HONGRAY GROUP

South Tongda Rd., East Dist. Jinzhou City Hebei, 052260, China

PRODUCT SPECIFICATION SHEET

STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

(Every mentioned standard is used in the latest edition)

Certification:	ISO 9001, ISO13485	
Conformity to guidelines:	Medical Device Regulation (EU) 2017/745: Class I under Rule 1 and Rule 5 according to Annex VIII (Subject to EU Declaration of Conformity)	
	Personal Protective Equipment Regulation (EU) 2016/425: Category III (Subject to EU Type-Examination Certificate)	
Conformity to standards:	Medical	PPE
	EN455	CAT III
	ASTM D6319	EN ISO21420
		EN ISO374
		EN 16523-1
Conformity to food regulation:	AP (2004) 4, (EU) No 10/2011	

DOCUMENT NO.: SETINO-NG-002

INITIALED DATE: 2021.10.25

REVISION DATE: N/A

REVISION TIMES: 0